



Five Star
ORTHODONTIC

Laboratory & Supply

P.O. Box 888

Denton, TX

76202-0880

1-800-521-2351

www.fivestarortho.com

EMA[®] - Custom

Simple ♦ Effective ♦ Comfortable



The EMA[®] - Custom appliance is a simple, patient-friendly oral appliance created for noninvasive treatment of snoring and OSA. The primary treatment mechanism of opening the bite and gently moving the mandible forward is achieved with the use of interchangeable elastic straps that offer varying degrees of mandibular advancement. The flexibility of these elastic straps provides unsurpassed lateral movement and overall TMJ comfort. The 2 mm thick pressure formed bases offer orthodontic retention (**resulting in no tooth movement**) and maximum anterior tongue space because there are no projections in the palate. The EMA[®] - Custom is available at [Five Star Orthodontic Lab & Supply](#) and has received FDA approval for marketing.

STRAPS		DUROMETER (FIRMNESS)			
		Soft (White)	Medium (Yellow)	Firm (Blue)	Extra Firm (Clear)
L E N G T H	#1 Extra Long 21 mm (Advances Mandible 4 mm)				
	#2 Long 19 mm (Advances Mandible 6 mm)				
	#3 Medium 17 mm (Advances Mandible 8 mm)				
	#4 Short 15.5mm (Advances Mandible 9.5 mm)				
	#5 Extra Short 14 mm (Advances Mandible 11 mm)				

Lab places buttons 25 mm apart